

FUJIIRYŌKI

Made in Japan



CYBER-RELAX

THE PIONEER IN MASSAGE CHAIRS

Since 1954

IT, ergonomics and craftsmanship coalesce to provide a massage that will blow you away.


Delivers the optimal massage for the sitter

Unique Function

3D point navigation system


The height, weight and shape of the human body as well as its lines change with time. For instance, there is a difference in the line of the back just after waking up in the morning and just after coming home from work in the evening.

The massage chair by Fujiiryoki is equipped with a unique function called the 3D Point Navigation Sensor. It automatically reads the sitter's shape of the back and begins massaging after determining the pressure points that are individualized for the sitter.



Unique Function

極メカ
きわみ
KIWAMIMECHA



Relieves tightness.

Fujiiryoki's massaging ball mechanism called KIWAMIMECHA (Super-Mechanism) uses two balls, which provide a powerful massaging action.

Its purpose-built processor makes it possible to deliver movement that is more like manual massage.

By finely moderating the strength of the massage and moving stereoscopically, it provides deep and powerful stimulation to ease away tightness from the neck to the lower back.


Neck massage

The massage balls improve blood circulation around the neck by gripping and stretching the nape.



Super Knead / Super Tapping

The tightness in the trapezius is relieved by providing deep stimulation to the upper part of the shoulders.



Super Back Knead / Super Back Tapping / Buttocks Massage

These functions are effective in soothing the soreness caused by the strain in the back or buttocks, neuralgia and muscle ache.



Stimulation provided by gentle air envelopes the whole body

Air Massage

Airbags placed throughout the chair for providing massage to the whole body from the shoulders to hands and arms, the back, calves and soles, improve the circulation of blood by delivering a variety of massages such as compression, squeeze and stretch.

You will feel pleasure from the gentle enveloping stimulation, the gripping and stretching action, and experience a sense of abandonment when the squeezed blood vessels are released - Enjoy the comfort made possible by the crystallization of Fujiiryoki's technologies.

The satisfying full-specification model offers graceful quality and ultimate comfort.

EC-3700

Unique Function

Making a diverse and sophisticated air-massage possible

エアマジック Air Magic

The high-capacity, high-output compressor effortlessly controls 41 air bags to provide Fujiiryoki's unique and innovative massage. The whole body is soothingly massaged by complex actions such as stretch and twist.



Combination of mechanism and air **734** types of massage

The combination of Fujiiryoki's Super-Mechanism with two massage balls and air massage makes it possible to provide 734 types of fine massage from the nape to the soles.

Unique Function

KIWAMIMECHA Massage




Neck massage

The massage balls improve blood circulation around the neck by gripping and stretching the nape.



Super Knead / Super Tapping

The tightness in the trapezius is relieved by providing deep stimulation to the upper part of the shoulders.



Super Back Knead / Super Back Tapping / Buttocks Massage

These functions are effective in soothing the soreness caused by the strain in the back or buttocks, neuralgia and muscle ache.



The whole body free of constraints thanks to wonders of 'Air Magic'. A feeling of relief for both body and soul.

NEW Total Body Stretch Massage

The airbags around the shoulders and calves work in concert to provide a pleasant and soothing stretching massage.

The shoulder airbags hold the shoulders.



The ottoman airbags grip the ankles and calves, stretching and pulling down.

The lower back is sterically stretched by the cross & twist action


Twist Massage

The air bags located in the area of the chair for the lower back and buttocks independently move from side to side, up and down to provide a 3D massage. It loosens tense muscles in the lower back and buttocks.



● Air Cross

Air to the lower back and buttocks is circulated diagonally to provide a firm stretching action.



● Air Twist

Air to the lower back and buttocks on one side of the body is circulated simultaneously to provide a gentle stretching action.

Shoulder Tip Massage

The shoulder tips where tension tends to build and symptoms of stiffness and fatigue appear are softly and gently massaged by the high-capacity double air bags to soothe away the stiffness and fatigue in the neck and shoulders.



Hand & Arm Massage

There is a concentration of pressure points in the hands and arms. Tension in these points is properly released by two types of massage: Grip and Knead.




Buttocks Air Massage

Soft air stimulation gently loosens the muscles in the buttocks and the back of the thighs to promote a better blood circulation in the groins.



Lower Back & Buttocks Massage

The interlocking action of intensive massage to the lower back and air massage to the buttocks makes it possible to provide relief to the lowermost part of the sidebone, which cannot be reached by ordinary massage.



Leg Massage

The air bags aid the "milking action", which promotes the circulation of blood in the lower body, by repeating compression and release from the lower part of the calves to the upper part of the calves, much like milking a cow.



Sole Massage

The soles are stimulated by air bags imbedded with pressure balls as the feet are firmly held from the heels.





Remote Control

Sophisticated functions can be operated simply by using the remote control with easy-to-understand user-friendly operating buttons.

Easy Quick Start

Quick Start a session by choosing your height from 3 different settings, or 'My Memory' in which you can register your course.

Register your course in 'My Memory'

Register a course individualized to you according to the shape of your back and the positions of your shoulders and lower back.

Cross-shaped Key for a simple operation

Choices can be easily made from a variety of functions by the use of a user-friendly cross-shaped key.

CYBER-RELAX



EC-3700 BK(ブラック)

CS(ベージュ) EC-3700



Electric Reclining

Find a comfortable, desired position via the easy to use remote control.



Perfect Fit Footrest Unique Function

The new length adjustable footrest, via a 2-step slide mechanism, ensures a massage in a comfortable position even for people of large build.



Accessories

Removable Back Pad & Pillow

The pad and pillow can be removed for easy cleaning.



Specifications

- Measurements (Approx): [When ottoman is stored] W750×D1,240×H1,210(mm) [In reclining position] W750×D2,010×H780(mm)
 - Reclining angle (Approx): 120 ~170 degrees
 - Weight (Approx): 86kg
 - Voltage: AC110V,120V,220V,230V,240V(50/60Hz)
 - Power consumption: 115W
 - Rated time: 30 minutes
- EC-3700's Japanese model AS-830 is 220AGBZX00122000 certified for medical device (Managed-Care Device) in Japan.

Model Number / Exterior	Upholstery	Remote Control	3D Point Navigation System	3D Mode Massage	Extreme Massage	Extreme Hip Massage Alignment	Fluctuation Massage	Waist/Hip Massage	Twist Massage	Sliding Footrest	Perfect Fit Footrest	Lower Body Stretch		Air Massage			Sole Pressure Balls	Auto-course	Massage Types	Mechta Massage Speed Adjustment	Width Adjustment	Air Massage Adjustment	Memory Function	Skip Function	Repeat Function	Pulse Mode	Rhythm Mode	Simultaneous Air Massage	Remote Control Stand	Airbags	Reclining	Dimensions (mm) (Reclined Position)	Weight	Power Consumption
												Foot Stretch	Double Stretch	Shoulder	Arm	Back / Waist																		
EC-3700	Leather (Black / Beige)	 ● User-friendly display ● Guidance Function ● Cross Key	●	●	●	●	●	●	●	●	●	●	●	●	● (new method)	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	W750×D1240×H1210 <W750×D2010×H780>	86 kg	115 W
EC-2700	Leather (Black)	 ● User-friendly display ● Guidance Function ● Cross Key	●	●	●	●	●	●	●	●	●	●	●	●	● (new method)	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	W750×D1300×H1170 <W750×D2040×H720>	81 kg	115 W
EC-1700	Leather (Black)	 ● Guidance Function ● Cross Key	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	W710×D1180×H1200 <W710×D1780×H720>	79 kg	115 W	

The standard model in a stylish design
equipped with advanced functions

EC-2700

Unique Function

Making a diverse and sophisticated
air-massage possible

エアマジック Air Magic

The high-capacity, high-output
compressor effortlessly controls
35 air bags to provide Fujiiryoki's
unique and innovative massage.
The whole body is soothingly
massaged by complex actions
such as stretch and twist.



Combination of
mechanism and air

677 types of massage

The combination of Fujiiryoki's Super-Mechanism with two massage balls and air massage makes it possible to provide 677 types of fine massages from the nape to the soles.

Unique Function

KIWAMIMECHA Massage



Neck massage

The massage balls
improve blood circulation
around the neck by
gripping and stretching
the nape.



Super Knead / Super Tapping

The tightness in the
trapezius is relieved
by providing deep
stimulation to the upper
part of the shoulders.



Super Back Knead / Super Back Tapping / Buttocks Massage

These functions are
effective in soothing
the soreness caused
by the strain in the back
or buttocks, neuralgia
and muscle ache.



Firm grip and slow stretch

Stretch Massage

The calves are massaged by the air bags in the ottoman by gripping and slowly stretching them. This massage improves the circulation of the blood in the lower body, which is prone to swelling.



● Stretch

The legs are held firmly and stretched forward.



● Leg Stretch

The legs are held firmly and stretched downward.

● Double Leg Stretch

The legs are held firmly, stretched forward and then downward.

The lower back is sterically stretched by the cross & twist action

Twist Massage

The air bags located in the area of the chair for the lower back and buttocks independently move from side to side, up and down to provide a 3D massage. It loosens tense muscles in the lower back and buttocks.



● Air Cross

Air to the lower back and buttocks is circulated diagonally to provide a firm stretching action.



● Air Twist

Air to the lower back and buttocks on one side of the body is circulated simultaneously to provide a gentle stretching action.

Hand & Arm Massage

There is a concentration of pressure points in the hands and arms. Tension in these points is properly released by two types of massage: Grip and Knead.



Leg Massage

The air bags aid the "milking action", which promotes the circulation of blood in the lower body, by repeating compression and release from the lower part of the calves to the upper part of the calves, much like milking a cow.



Buttocks Air Massage

Soft air stimulation gently loosens the muscles in the buttocks and the back of the thighs to promote a better blood circulation in the groins.



Sole Massage

The soles are stimulated by air bags imbedded with pressure balls as the feet are firmly held from the heels.



Lower Back & Buttocks Massage

The interlocking action of intensive massage to the lower back and air massage to the buttocks makes it possible to provide relief to the lowermost part of the sidebone, which cannot be reached by ordinary massage.



Bringing peace and comfort to life.
The enhanced-function easy-going model.

EC-1700



Combination of mechanism and air **338** types of massage

The basic actions of the massage balls combined with air massage have made it possible to deliver a more varied and fine 338 types of massage.

Unique Function

KIWAMIMECHA Massage




Neck massage

The massage balls improve blood circulation around the neck by gripping and stretching the nape.



Super Knead / Super Tapping

The tightness in the trapezius is relieved by providing deep stimulation to the upper part of the shoulders.



Super Back Knead / Super Back Tapping / Buttocks Massage

These functions are effective in soothing the soreness caused by the strain in the back or buttocks, neuralgia and muscle ache.



Unique Function

Loosening Massage





Knead & Tapping Automatic Variable Function

A new sensation a machine could not provide before. A manual-like massage makes it possible to provide a fine and rhythmical massage. Tension is gently released by the changeable rhythm instead of a steady mechanical rhythm.



Knead
The kneading sensation feels very much like a real massage.



Tapping
By controlling the speed of pounding and changing it gradually, it provides an effective massage as the body acclimatizes to it.

Massage at your convenience without having to take your shoes off.

Leg Massager

It promotes a better circulation of blood in the legs by firmly taking hold of the calves with the air bags. It is equipped with a variety of functions such as its unique Cross-press system.



Unique Function

Milking Action

The milking action sends blood back to the heart by pressurizing the blood vessels. Lack of exercise and staying in the same position for many hours cause the muscles in the legs to weaken, resulting in swelling or tiredness.

Legs are firmly held and slowly stretched

Leg Stretch

The Stretch Massage uses the air bags to hold the calves and slowly pulling them downward. It promotes a better circulation of blood in the legs, which are prone to swelling.





Remote Control

A user-friendly remote control made from a user's standpoint.

Large letters and icons are displayed clearly on the 3.4-inch LCD display.

Anyone can read and check the usage status easily for the icons are easy to recognize and the letters are displayed large.

Cross-shaped Key for a simple operation

Choices can be easily made from a variety of functions by the use of a user-friendly cross-shaped key.

Operate the variety of functions without confusion

The key layout has been centralized according to the types of function so that the user can operate the variety of advanced functions without confusion.

CYBER-RELAX



EC-1700 BK (black)

Electric Reclining

Find a comfortable, desired position via the easy to use remote control.



Specifications

●Measurement (Approx): [When ottoman is stored] W710×D1,180×H1,200(mm) [In reclining position] W710×D1,780×H720(mm) ●Reclining angle (Approx):120 ~170 degrees ●Weight (Approx): 79kg ●Voltage: AC110V,120V,220V,230V,240V(50/60Hz) ●Power consumption: 115W ●Rated time: 30 minutes EC-1700's Japanese model AG-660 is 220AGBZX00197000 certified for medical device (Managed-Care Device) in Japan.

Model Number / Exterior	Upholstery	Remote Control	3D Point Navigation System	3D Mode Massage	Extreme Massage	Extreme Hip Massage Alignment	Fluctuation Massage	Waist/Hip Massage	Twist Massage	Sliding Footrest	Perfect Fit Footrest	Double Stretch	Lower Body Stretch		Air Massage		Sole Pressure balls	Foot Stretch	Sole Stretch	Massage Types	Auto-course	Quick-course	Mecha Massage Adjustment	Mecha Massage Speed Adjustment	Width Adjustment	Air Massage Adjustment	Memory Function	Skip Function	Repeat Function	Pulse Mode	Rhythm Mode	Simultaneous Air Massage	Remote Control Stand	Airbags	Reclining	Dimensions (mm) (Reclined Position)	Weight	Power Consumption		
													Shoulder	Arm	Hip / Thigh	Back / Waist																							Shoulder	Arm
EC-3700	Leather (Black / Beige)	 ● User-friendly display ● Guidance Function ● Cross Key	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	7 3 4	1 2	3	7 Step (7 Step Adjustment in Course Mode)	3 0 ~ 3 3	3 0 ~ 3 3	3 0 ~ 3 3	5 Step (5 Step Adjustment for Arm)	●	●	●	●	●	●	●	●	4 1	120-170° (Approximate) HydroPump Chair back and footrest Independently adjustable footrest (up/down)	W750×D1240×H1210 <W750×D2010×H780>	8 6 kg	1 1 5 W
EC-2700	Leather (Black)	 ● User-friendly display ● Guidance Function ● Cross Key	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	6 7 7	1 2	3	7 Step (7 Step Adjustment in Course Mode)	3 0 ~ 3 3	3 0 ~ 3 3	3 0 ~ 3 3	5 Step (5 Step Adjustment for Arm)	●	●	●	●	●	●	●	●	3 5	120-170° (Approximate) HydroPump Chair back and footrest Independently adjustable footrest (up/down)	W750×D1300×H1170 <W750×D2040×H720>	8 1 kg	1 1 5 W
EC-1700	Leather (Black)	 ● Guidance Function ● Cross Key	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	3 3 3	1 0	1	7 Step (5 Step Adjustment in Course Mode)	3 0 ~ 3 3	3 0 ~ 3 3	3 0 ~ 3 3	5 0 ~ 3 3	●	●	●	●	●	●	●	●	1 0	120-170° (Approximate) HydroPump Chair back and footrest Independently adjustable footrest (up/down)	W710×D1180×H1200 <W710×D1780×H720>	7 9 kg	1 1 5 W

Half a century of making massage chairs FUJIIRYOKI is an all-round manufacturer for health and beauty

History of massage chairs

1954

The first FUJIIRYOKI massage chair

The world's first wooden massage chair. The founder, Nobuo Fujimoto, commercialized it after creating numerous prototypes made from scraps and junk. It still works when it is plugged in.



1965

FUJIIRYOKI Automatic Massage Chair L-11

The chair uses a 4-ball system, which can provide an extensive massage. It became a hot-seller as the popular version of massage chairs. The height can be adjusted by the handles on the sides and the strength of the massage can be modulated by a switch.



1975

Massage Chair FM-52

Fujiiryoki's unique special clutch made it possible to provide innovative fully automated 'multiple massage functions' such as kneading and pounding.



1987

Roller Chair KIT-10DX

The rolling massage chair used 8 rollers. It was a hot seller for it had significantly improved features such as remote control and spurred rollers.



1995

Air Massage Chair Royal Chair MC-133

The industry's first massage chair with air bags. It was especially innovative for its massage functions for the lower body (thighs and legs), which were not available with conventional massage chairs.



2000

Massage Chair Air Solution AS-777

The 'air bags' and 'massage balls' evolved in their own ways. It was a new idea to combine two types of massage. This chair was a concept model, which determined the trend of massage chairs in the 21st century.



2001

Massage Chair CYBER-Relax AS-001

This model was equipped with the industry's first, many groundbreaking functions such as the 3D Point Navigation System, which senses the body shape and automatically chooses the best massage program for the sitter.



2002

Massage Chair CYBER-Relax AS-003

This model was created with organic and natural flow in mind, and was also highly regarded for its design. It was particularly well received for its air massage functions for the legs, including 'Leg Stretch' and 'Cross Press.'



2004

Massage Chair CYBER-RELAX AS-800

The first massage chair to be equipped with functions for the shoulder tips and soles. It made it possible to provide an extensive massage from the nape to the soles.



2008

Massage Chair CYBER-RELAX AS-830

The latest flagship model is filled with all the details Fujiiryoki cares about. The combination of mechanical and air massages makes it possible to provide 734 types of massage such as the 'whole body stretch.'



About FUJIIRYOKI

Established in 1954, Fujiiryoki is a manufacturer of medical device, which provides health and beauty-related products. We boast the largest share (over 30%) in the high value-added massage chair market, and we have been distributing products to more than 20 countries throughout the world. In recent years, we recorded high-growth in the training equipment sector. We have also been supporting the health and beauty of customers worldwide for many years through our medical device such as potential, ultrashort wave and low-frequency therapy equipment and hearing aids.

Corporate History

1954	<ul style="list-style-type: none"> ● Fuji Seisakusho is established in Hannancho, Abeno-ku, Osaka City, to produce healthcare equipment. ● The first wooden massage chair is developed.
1960	<ul style="list-style-type: none"> ● The industry's first 4-ball massage chair is developed.
1965	<ul style="list-style-type: none"> ● The company reorganizes as a corporate body under the name Kabushikigaisha Fujiiryoki. ● The head office and factory move to Sugimotocho, Sumiyoshi-ku, Osaka City. ● Fujiiryoki Massage Chair's massage method "pounding" is patented in the USA, and obtains exclusivity rights in the US.
1971	<ul style="list-style-type: none"> ● Sakai factory is built and put into operation.
1982	<ul style="list-style-type: none"> ● The new office building is completed in Nipponbashi, Naniwa-ku, Osaka City.
2003	<ul style="list-style-type: none"> ● CYBER-Relax S.O massage chair series wins the Good Design Special Prize.
2004	<ul style="list-style-type: none"> ● CYBER-Relax AS-800 massage chair series wins the Good Design Award.
2005	<ul style="list-style-type: none"> ● The company obtains the ISO13485:2003 certification, which is the international standard for 'quality management system for healthcare equipment.'
2008	<ul style="list-style-type: none"> ● The head office moves to Chuo-ku, Osaka City. ● The Tokyo Showroom is opened.
2009	<ul style="list-style-type: none"> ● The Osaka factory and new office building are completed in Taishicho, Minamikawachi-gun, Osaka Prefecture.



Safety
Precautions

- Please read the [User Manual] before use and use this product properly.
- The following persons are strongly recommended not to use this product. Certain physical abnormalities may occur.
 - Persons prohibited from receiving a massage by a medical physician, for example those suffering from thrombosis, severe aneurysm, acute varicose veins, various kinds of dermatitis and skin infections (including inflammations of subcutaneous tissue), etc.
- The following persons must consult a medical physician before using this product as certain physical abnormalities may occur.
 - Persons with malignant tumors. ● Persons who are pregnant or menstruant, persons thought to be pregnant, or those immediately after childbirth.
 - Persons who suffer from paresthesia related to peripheral circulatory disturbance due to diabetes etc. ● Persons with body temperature exceeding 38°C (when strong acute inflammation symptoms, malaise, chills, blood pressure fluctuations are observed, when debilitated) ● Persons with wounds in treatment areas ● Persons who suffer from acute illness (dolorific) ● Persons who suffer from osteoporosis ● Persons who suffer from spinal column fractures ● Persons suffering from inflammation such as sprains and muscle strains ● Persons who suffer from paresthesia ● Persons who suffer from heart trouble ● Persons using a medical electronic device embedded in the body such as a pacemaker susceptible to electromagnetic disturbances ● Persons with abnormalities in the spine ● Persons who suffer from anemia ● Persons who need to massage a part of the body treated medically in the past ● Persons who need bed rest ● Persons receiving medical treatment ● Persons who feel an abnormality on the body other than the above
- Please do not allow unattended children or persons unable to declare their own intentions to use this chair. Please do not allow infants near the chair.

● Concerning the use of this catalogue

- The material contained in this catalogue are based on domestic Japanese Pharmaceutical Affairs Laws. The claims made in this catalogue may not be valid or relevant in your country.
- Fujiiryoki will not be held responsible for the use of this catalogue outside of Japan, and all responsibility will lie with the user in such circumstances.
- Product images shown may differ from actual product appearance.
- Designs and specifications are subject to change without notice.

FUJIIRYŌKI

Fuji Medical Instruments Mfg. Co., Ltd.
Overseas Sales Department

14th Floor, Oe Building, 1-22, Noninbashi 1-chome, Chuo-ku, Osaka, 540-0011, Japan
Tel: +81-6-4793-0745 Fax: +81-6-4793-0781



RoHS



FDA

